



CORPORATE PRESENTATION 2008

Member of the Lynk Biotech Group





Our Services

We specialize in the contract manufacturing of semi-solid topical dosage forms (creams, lotions, balms, gels and ointments).

Other services:

Tube-Filling

Secondary Packaging

Procurement Warehousing

Quality Testing

Warehousing

Inventory Management





Our Facility

- **Approx. 1,000 sqm manufacturing facility**
- **Compliant with PIC/S Pharma Good Manufacturing Practice (GMP) standards**
- **Fully equipped with a Class 100,000 clean room (PIC/S Grade D)**



We're Well Equipped

- Pilot 50-kg mixer
- Fully automated 500-kg mixer
- Ink jet coding machine
- Tube filling machine

All mixers are double jacketed and come with vacuum and homogenization as options.

Half-Ton Romaco Mixer



- Equipped with heater and chiller to enable hot and cold processes.
- Fast and efficient cleaning with built-in CIP/WIP (Clean-In-Place / Wash-In-Place).
- Efficient heat exchange for minimum processing time.
- Strong vacuum for optimal de-aeration.

50-kg Pilot Mixer



- Manufacturing vessel comes with a jacket that is suitable for water heating and water cooling at 2 kg/cm²
- Capable of performing CIP with the presence of 2 spray balls and a high pressure pump.
- Built-in Homogenizer with a maximum speed of 1440 rpm. Strong vacuum for optimal de-aeration.
- Built-in anchor agitator with a speed range of 10 to 50 rpm.

● ● ● Auto Tube Filling Machine



- Tube Diameter Range
- 10 to 52 mm
- Tube Length Range
- 50 to 250 mm
- Dosing Range
- 1.5 to 400 ml
- Output Range
- up to 70 tubes/min

EP Water Purification System



- Total viable aerobic count
- ≤ 100 cfu/ml
- Total organic carbon limit
- ≤ 0.5 mg/l
- Heavy metals
- ≤ 0.1 ppm
- Nitrates
- ≤ 0.2 ppm
- Conductivity at 25°C
- ≤ 5.1 $\mu\text{S/cm}$

Class 100K Dispensing Booth



- Top Filter Unit containing two pieces of HEPA filters with an efficiency of 99.99% @ 0.3 μm .
- Middle Filter Unit containing secondary filters; 60-65% efficiency.
- Bottom Filter Unit containing pre-filters; 25-30% efficiency.

Contact Details

Lynk Biotechnologies (Head Office)

Greatland Building
48 Lorong 21 Geylang #03-01
Singapore 388464

Tel: +65 6745 1550
Fax: +65 6745 5475
Email:
custsvc@lynkbiotech.com

www.lynkbiotech.com

Dynalynk Pharma

23 Gul Avenue
Singapore 629663

Tel: +65 6862 9168
Fax: +65 6862 9160
Email:
enquiry@dynalynkpharma.com

www.dynalynkpharma.com



Thank You



**Your Trusted & Reliable
Manufacturing Partner**